

Capre F Treatment table

Versatility for physiotherapists

The modifiability, adjustability and silent operation of the sturdy Lojer Capre F improve functionality and enhance the comfort of use. Its multiple height settings, easy-to-use armrests and side supports, the stepless adjustable extension and flexion angle of the foot section improve the ergonomics and efficiency of the work being carried out. All Lojer Capre treatment tables come with high-quality, durable FUTURE upholstery.

*Lojer Capre – Leading
the Way in Treatment
Table Safety and Design*



**Lojer Capre F3
with optional accessories**



**Lojer Capre F5
with optional accessories**




Check out the new Capre product
video: <http://bit.ly/capre2nd>

LOJER® For easy care

Lojer Capre F	Capre F3	Capre F3H	Capre F5	Capre F5H
Table sections	3	3	5	5
Width (cm)	55,65,75	55,65,75	55,65,75	55,65,75
Upholstery FUTURE (4 options)	•	•	•	•
Upholstery FUTURE (10 additional options)	option	option	option	option
Anatomic design upholstery	option	option	option	option
Fully adjustable concave arm rests (gas springs)	-	-	•	•
Turning arm rest (can be turned to side supports)	option	option	option	option
Side supports with fast locking and horizontal adjustment	option	option	option	option
Head support adjustment with gas spring	•	•	•	•
Head section, oval shape	option	option	option	option
Head section, U-shape	option	option	option	option
Double hinged head support	-	-	option	option
Leg raise with gas spring	+70°/-12°	+70°/-12°	+70°/-12°	+70°/-12°
Drainage with gas spring	option	option	option	option
Drainage with electric adjustment	option	-	option	-
Central locking castors (Ø75mm)	option	option	option	option
Foot adjustment (levelling)	•	•	•	•
Height adjustment range (cm)	45-92	49-90	45-92	49-90
Electric height adjustment with foot pedal	•	-	•	-
Hydraulic height adjustment with foot pedal	-	•	-	•
Additional hand control or foot pedal	option	-	option	-
Height adjustment with 360° foot bar	option	-	option	-
Paper roll holder (50 or 60 cm)	option	option	option	option
Battery operation	option	-	option	-
Safe working load (SWL)	210 kg	210 kg	210 kg	210 kg
CE-marking according to 93/42/EEC	•	•	•	•
Warranty (years)	2/10 y	2/10 y	2/10 y	2/10 y
Country of Manufacture	Finland	Finland	Finland	Finland

Explanation of markings: • =standard feature option =optional accessory - =not available



Capre F5 fully adjustable arm rests



Capre F3/F5 the stepless adjustable extension and flexion angle of the footsection



Drainage with gas spring or electric motor



Side supports with fast locking and horizontal adjustment



Check out the new feature videos of Capre:
<http://bit.ly/caprefeatures>



Upholstery options F3, F5, F5 hourglass



LOJER® For easy care

This document template PC001 has been approved 01.01.2016

MADER
 MEDICAL SYSTEMS

Kirchlergründe 6/19
 A-6200 Jenbach

T: +43 5244 211 41 www.mader-medical.com
 F: +43 5244 211 41-42 office@mader-medical.com